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










Made in France  
 ACCU-CHEK, ACCU-CHEK SPIRIT COMBO and DIAPORT are trademarks of Roche.

© 2020 Roche Diabetes Care

 Roche Diabetes Care GmbH  
 Sandhofer Strasse 116  
 68305 Mannheim, Germany  
 www.accu-chek.com

## EN Accu-Chek DiaPort fixation disc

### Intended Use

The Accu-Chek DiaPort system is intended for continuous intraperitoneal insulin infusion (CIPI). For this purpose, it shall be used together with an external Accu-Chek Spirit Combo Insulin pump. The Accu-Chek DiaPort system delivers insulin from the Accu-Chek insulin pump through the infusion set and the catheter into the peritoneal cavity.

### ⚠ WARNING

#### Risk of hyperglycaemia

The use of a non-compatible insulin pump can lead to insufficient insulin delivery. Use the Accu-Chek DiaPort system only with an Accu-Chek Spirit Combo insulin pump.

Intended users of the Accu-Chek DiaPort system are insulin-dependent people with type I or type II diabetes who meet the indication for CIPI therapy. The Accu-Chek DiaPort system is only for use by patients and healthcare professionals specifically trained in the use of the individual components relative to their areas of expertise. The Accu-Chek DiaPort fixation disc helps to keep the port body in place.

### Indications and Contraindications

See the Accu-Chek DiaPort Handbook for Port Users.

### Component Overview

This is a sterile, single use product that has been sterilised with ethylene oxide. The fixation disc and the disc applicator must be replaced at least every 6 days, for example, when changing the cartridge of the insulin pump.

### ⚠ PRECAUTION

#### Risk of infection

The use of system components beyond the specified replacement intervals can cause infections. Replace system components according to the specified intervals.

- A Fixation disc
- B Disc applicator

### ⚠ WARNING

#### Risk of infection

A damaged packaging may no longer be sterile or may contain a damaged component. Before using a new component, make sure its packaging is not damaged. If the packaging is damaged, do not use its contents.

#### Risk of different types of harm

Modification of system components or disregard of the instructions for use may prevent the system from working as intended. This can cause different types of harm, for example, infection, hyperglycaemia, damage to the skin, or scarring due to explantation of the system.

Do not modify the system and adhere to the instructions for use. **Risk of suffocation** This product contains small parts. Keep small parts away from people who might swallow them, for example, small children.

### Storage

Keep the product dry and away from sunlight. Store the product between 5 and 45 °C. Store the disc applicator in a clean and dry environment.

### ⚠ WARNING

#### Risk of hyperglycaemia and infection

The use of sterile components after their *Use by* date can lead to infections and to abscesses. Do not use components that have passed their *Use by* date.

### Attaching the Fixation Disc

- 1 Wash and dry your hands thoroughly. Clean the port with sterile sodium hypochlorite solution or clean water. Do not use solutions containing alcohol, as these can be too aggressive for the skin.
- 2 Close the front arms of the disc applicator to attach the disc applicator to the fixation disc. The connection is correct when there is no gap between the disc applicator and the fixation disc.
- 3 Close the rear arms of the disc applicator to spread the fixation disc.
- 4 Position the fixation disc on the port. Turn it until it locks. The placement is correct when the opening of the fixation disc is no greater than 1 mm (0.04 inches).
- 5 Close the front arms of the disc applicator to remove the disc applicator.

### ⚠ PRECAUTION

#### Risk of hyperglycaemia

Handling errors or defective components can lead to insufficient insulin delivery. After replacing a component of the system, check your blood glucose level at least once within 1 to 3 hours.

### Risk of adverse reaction

In rare situations, the implanted system may cause adverse reactions. If you experience an adverse reaction, contact your healthcare professional immediately.

### Removing the Fixation Disc

- I Wash and dry your hands thoroughly. Use the handling aid to disconnect the infusion set from the port membrane.
- II Close the front arms of the disc applicator to attach the disc applicator to the fixation disc.
- III Close the rear arms of the disc applicator to spread the fixation disc.
- IV Remove the fixation disc.

### ⚠ WARNING

#### Risk of interruption of therapy and scarring

The fixation disc prevents the port body from sinking into the subcutaneous tissue, which may require explantation. Wear the fixation disc at all times.

#### Risk of interruption of therapy and scarring

If the port body sinks into the subcutaneous tissue, explantation may be required. If you notice that the port body starts sinking, contact your healthcare professional immediately.

### Cleaning

Wash and dry the fixation disc daily. Do not place cream around the port area.

### Disposal

Discard the used product according to local regulations. For more information, contact your local council or authority.

### Reporting of Serious Incidents

For a patient/user/third party in the European Union and in countries with identical regulatory regime; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

### Summary of Safety and Clinical Performance

You find the SSCP after the launch of the European database on medical devices (Eudamed) here: <https://ec.europa.eu/tools/eudamed>

### Distributed in the United Kingdom by:

Roche Diabetes Care Limited  
Charles Avenue, Burgess Hill  
West Sussex, RH15 9RY, United Kingdom  
Accu-Chek Pump Careline 18  
UK Freephone number: 0800 731 22 91  
ROI Freephone number: 1 800 88 23 51

### 1) calls may be recorded for training purposes

Some mobile operators may charge for calls to these numbers.

[burgesshill.insulinpumps@roche.com](mailto:burgesshill.insulinpumps@roche.com)  
[www.accu-chek.co.uk](http://www.accu-chek.co.uk)  
[www.accu-chek.ie](http://www.accu-chek.ie)

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Roche Diabetes Care France SAS