Handbook for Diabetologists, Surgeons and Diabetes Nurses

The Accu-Chek DiaPort system

EN DE FR IT NL NO

Ⓔ	Handbook for Diabetologists, Surgeons and Diabetes Nurses	2-39
Ⓓ	Handbuch für Diabetologen, Chirurgen und Diabetesberater	40-77
Ⓕ	Guide pour diabétologues, chirurgiens et infirmiers en diabétologie	78-115
Ⓘ	Manuale per diabetologi, chirurghi e infermieri specializzati nell'assistenza dei pazienti con diabete	116-153
Ⓝ	Handboek voor diabetologen, chirurgen en diabetesverpleegkundigen	154-191
Ⓝ	Håndbok for diabetologer, kirurger og diabetessykepleiere	192-229

Distributed in the United Kingdom by:

Roche Diabetes Care Limited
Charles Avenue, Burgess Hill
West Sussex, RH15 9RY, **United Kingdom**
Accu-Chek Pump Careline ¹⁾:
UK Freephone number: 0800 731 22 91
ROI Freephone number: 1 800 88 23 51

¹⁾ calls may be recorded for training purposes
Some mobile operators may charge for calls to these numbers.

burgesshill.insulinpumps@roche.com
www.accu-chek.co.uk
www.accu-chek.ie

Australia

Roche Diabetes Care Australia Pty. Limited
Insulin Pump Support: 1800 633 457
australia.insulinpumps@roche.com
www.accu-chek.com.au

Product Overview

Components of the Accu-Chek DiaPort system:

<p>1</p> 	<p>Handling aid</p>
<p>2</p> 	<p>Infusion set with ball cannula</p>
<p>3</p> 	<p>Disc applicator</p>
<p>4</p> 	<p>Fixation disc</p>

5

Membrane


Port body

6

Catheter

7

Implant card

About this Handbook

The purpose of this handbook is to support healthcare professionals and surgeons in the care of people with diabetes who consider the Accu-Chek DiaPort system as an option for insulin delivery. A separate handbook is available for patients. The Accu-Chek DiaPort is referred to as “port” when talking about the port body that contains the membrane and as “system” when talking about the port and its accessories, such as, the catheter or the fixation disc.

This handbook is intended to:

- Provide you with essential information about the therapeutic aspects of switching from continuous subcutaneous insulin infusion (CSII) to continuous intraperitoneal insulin infusion (CIPII).
- Help you to understand what is involved in the implantation of the system and post-operative care.
- Give you an insight into the patient’s daily life with the port after the operation.

The Accu-Chek DiaPort system is only for use by healthcare professionals and patients specifically trained in the use of the individual components relative to their areas of expertise. The system may only be implanted by a trained surgeon working with a certified Accu-Chek DiaPort center of excellence.

Prior to using this product or any of its components, consult the instructions for use that come with the product. For explanations of the medical terms used in this handbook, see the Glossary section.


Table of Contents

Product Overview	3
About this Handbook.....	5
Introduction	7
1 Deciding on the Accu-Chek DiaPort System.....	8
1.1 Indications.....	8
1.2 Intended Users	9
1.3 Contraindications	9
1.4 General Warnings.....	10
1.5 Advantages and Disadvantages	11
1.6 Clinical Experiences	13
2 System Components.....	15
2.1 Port Body.....	15
2.2 Membrane.....	16
2.3 Catheter.....	17
2.4 Fixation Disc	18
2.5 Infusion Set.....	19
3 Implantation.....	20
3.1 Implantation Set.....	20
3.2 Surgical Procedure	21
3.3 Implantation Site.....	22
3.4 What Will the Implanted Port Look Like?.....	23
4 Postoperative Care	24
5 Therapy.....	25
5.1 Insulin Kinetics With CIPII Therapy.....	25
5.2 Switching from CSII to CIPII.....	26
5.3 Hypoglycaemia and Hyperglycaemia.....	28
5.4 Catheter Replacement.....	29
5.5 Membrane Replacement.....	31
6 Daily Care	32
7 Medical Conditions not Directly Related to the Port.....	33
8 Daily Life	34
Glossary	37
Symbol Descriptions	39

Introduction

The Accu-Chek DiaPort system enables continuous intraperitoneal insulin infusion: That is, infusion of insulin into the peritoneal cavity with a compatible Accu-Chek insulin pump and a specially designed ball cannula infusion set. The port has a titanium body with a catheter that is placed in the abdomen. The top of the port rises above the surface of the skin by about 5 mm (0.19 inches), while a flower-shaped plate is surgically placed under the skin in a subcutaneous pocket, providing stability for the system.

The system has been designed for people with diabetes who cannot fully benefit from continuous subcutaneous insulin infusion (CSII).


The Accu-Chek DiaPort system with an Accu-Chek Spirit Combo insulin pump.

1 Deciding on the Accu-Chek DiaPort System

1.1 Indications

The Accu-Chek DiaPort system may be a suitable solution if one or more of the following indications exist:¹

- People with diabetes mellitus requiring insulin in whom subcutaneous insulin therapy fails as shown by not achieving glycaemic targets, or only achieving glycaemic targets at the expense of increased rates of severe hypoglycaemia in spite of the use of advanced technology, for instance in subcutaneous insulin resistance.
- People with diabetes mellitus requiring insulin in whom subcutaneous insulin is not feasible because of skin disorders, for instance in the presence of lipohypertrophy or -atrophy.

¹ Schade D.S., Duckworth W.C., *In Search of the Subcutaneous-Insulin-Resistance Syndrome*. N Engl J Med, 1986 Jul 17;315(3):147-53.
Soudan B., et al., *Extreme Subcutaneous Insulin Resistance: A Misunderstood Syndrome*. Diabetes Metab, 2003 Nov;29(5):539-46.
Riveline J.P., et al., *Subcutaneous Insulin Resistance Successfully Circumvented on Long Term by Peritoneal Insulin Delivery From an Implantable Pump in Four Diabetic Patients*. Diabetes Metab, 2005 Nov;31(5):496-8.
Baillot-Rudoni S., et al., *Implantable Pump Therapy Restores Metabolic Control and Quality of Life in Type 1 Diabetic Patients With Buschke's Nonsystemic Scleroderma*. Diabetes Care, 2006 Jul;29(7):1710.
Hilgard D., Laengler A., *Severe subcutaneous insulin infusion resistance in a 13 year old child - a case-report of successfully performed therapy with a diaport™-system*. Pediatric Diabetes, 2009; 10(Suppl 11):37.
Noud M.N., et al., *Benefits of intra-peritoneal insulin administration in a child with severe insulin-induced lipoatrophy*. Pediatric Diabetes, 2009;10 (Suppl 11):100.
Lee S.I., Narendran P., *Intraperitoneal Insulin Therapy for a Patient With Type 1 Diabetes With Insulin Injection Site Inflammation*. BMJ Case Reports, 2014 Aug 21.
Garcia-Verdugo R., Erbach M., Schnell O., *A New Optimized Percutaneous Access System for CIPII*. J Diabetes Sci Technol. 2017 Jul;11(4):814-821..

1.2 Intended Users

Intended users of the Accu-Chek DiaPort system are insulin-dependent people with type I or type II diabetes who meet the indication for CIPII therapy.

The user or his/her permanent caregiver must be fully aware of the details of the Accu-Chek insulin pump therapy and he or she must be well-trained in the use of the insulin pump therapy as well as the backup systems.

The user's healthcare professional decides whether a user, either alone or with mandatory and permanent support of a caregiver is able to use an insulin pump for intraperitoneal infusion. Special attention has to be paid when using intraperitoneal insulin pump therapy in children or neonates. There is no experience with intraperitoneal insulin pump therapy in neonates.

WARNING

Risk of inappropriate therapy

The Accu-Chek DiaPort system is designed for intended users only. Do not use the system with patients that do not meet the indications, or where one or more contraindications apply.

1.3 Contraindications

The Accu-Chek DiaPort system may not be suitable if one or more of the following contraindications apply:

- Patients not considered suitable for using an external insulin pump, or patients with contraindications for CSII:
 - Poor compliance with current treatment
 - Recurrent failure to attend clinics
 - Inability to follow a treatment plan
 - Unwillingness to calculate meal dosages or to perform the necessary number of blood glucose tests per day for this type of therapy
 - Evidence of psychiatric conditions
- Allergy to one or more of the port materials (titanium, polyethylene, etc.) including its infusion set
- Disorders of the gastrointestinal tract interfering severely with the intraperitoneal route of administration (for example, colonic disease)
- Continuous ambulatory peritoneal dialysis (CAPD)
- Concomitant drugs administered by intraperitoneal infusion (for example, chemotherapy)
- Severely impaired immune response
- High levels of anti-insulin antibodies

1.4 General Warnings

Do not use this product for any medication other than U100 short-acting regular human insulin or sterile physiological saline solution.

Replace the infusion set and the fixation disc according to the intervals specified in the respective instructions for use.

PRECAUTION

Risk of hyperglycaemia

The use of system components beyond the specified replacement intervals can cause occlusions or leaks that can lead to insufficient insulin delivery.

Replace system components according to the specified intervals.

Risk of hyperglycaemia

Handling errors or defective components can lead to insufficient insulin delivery.

After replacing a component of the system, check the blood glucose level at least once within 1 to 3 hours.

WARNING

Risk of infection

A damaged packaging may no longer be sterile or may contain a damaged component.

Before using a new component, make sure its packaging is not damaged. If the packaging is damaged, do not use its contents.

Risk of suffocation

This product contains small parts.

Keep small parts away from people who might swallow them, for example, small children.

Risk of hyperglycaemia and infection

The use of sterile components after their *Use by* date can lead to infections and to abscesses.

Do not use components that have passed their *Use by* date.

Risk of different types of harm

Modification of system components or disregard of the instructions for use may prevent the system from working as intended. This can cause different types of harm, for example, infection, hyperglycaemia, damage to the skin, or scarring due to explantation of the system.

Do not modify the system and adhere to the instructions for use.

1.5 Advantages and Disadvantages

Advantages

Therapy using the Accu-Chek DiaPort system can have the following advantages:

- In contrast to subcutaneous injections or subcutaneous pump therapy, the insulin flows through the catheter into the peritoneal cavity and reaches the liver more quickly. The insulin starts taking effect within a few minutes.
- Using this method may improve your diabetes control, especially reducing frequent hypoglycaemia, as CIPII more closely mimics the normal physiological secretion of insulin. CIPII can lead to an improved blood glucose profile.
- Compared to subcutaneous insulin infusion, HbA_{1c} values can be improved.²
- The frequency of unexpected and rapid drops in the blood glucose level can be reduced.
- You no longer need to change infusion sites every 2 to 3 days, as the insulin is infused directly through the port catheter.

Disadvantages


Therapy using the Accu-Chek DiaPort system can have the following disadvantages:

- There is a risk of infection of the tissue around the top of the port body, especially just after the implantation. However, taking the recommended precautions after implantation helps to reduce the risk of infections.
- There can be adhesions to the abdominal wall, which may result in pain, or the tip of the catheter may become encapsulated. In this case, it is sometimes necessary to replace the catheter.
- The catheter can become blocked, which is called occlusion.
- When there are serious or recurring infections, the port must be explanted. However, once the infection is cleared, a new port can be implanted at a new site.

² Noud M.N., et al., *Benefits of intra-peritoneal insulin administration in a child with severe insulin-induced lipodystrophy*. Pediatric Diabetes, 2009;10(Suppl 11):100.

Liebl A., et al., *Long-term clinical evaluation of the new Accu-Chek® diapor, a port system for continuous intraperitoneal insulin infusion: 24-month results*. Diabetes, 2014; 63: A241.

The amount of hypoglycaemic and severe hypoglycaemic events in CIPII is reduced compared to CSII.³ One study showed that the number of incidents of severe hypoglycaemia with CIPII was less than half the number compared to CSII.

³ Liebl A., et al., *A Reduction in Severe Hypoglycaemia in Type 1 Diabetes in a Randomized Crossover Study of Continuous Intraperitoneal Compared With Subcutaneous Insulin Infusion*. *Diabetes Obes Metab*, 2009 Nov;11(11):1001-8.

Haak T., et al., *New path in intraperitoneal insulin administration*. *Diabetes Forum*, 2016; (6): 11-15, Kirchheim-Verlag

Garcia-Verdugo R., et al., *A New Optimized Percutaneous Access System for CIPII*. *J Diabetes Sci Technol*, 2017 Jul;11(4):814-821.

Reix N., et al., *Glycemic Management of Diabetes by Insulin Therapy*. *Minerva Endocrinol*, 2013 Mar;38(1):29-46.

Van Dijk P.R., Logtenberg S.J., Groenier K.H., Gans R.O., Kleefstra N., Bilo H.J., *Continuous Intraperitoneal Insulin Infusion in Type 1 Diabetes: A 6-year Post-Trial Follow-Up*. *BMC Endocr Disord*, 2014 Apr 7;14:30.

Spaan N., Teplova A., Stam G., Spaan J., Lucas C., *Systematic Review: Continuous Intraperitoneal Insulin Infusion With Implantable Insulin Pumps for Diabetes Mellitus*. *Acta Diabetol*, 2014;51(3):339-51.

Gehr B., et al., *Continuous intraperitoneal insulin infusion shows a great benefit-risk ratio under long-term real-world use in a vulnerable population with diabetes*. 2019. EASD Poster.

1.6 Clinical Experiences

Over the last 20 years, implantable pumps were the most frequently clinically used systems for continuous intraperitoneal insulin infusion.⁴ The EVADIAC group, a network of collaborative centers in France, collected experience with implantable pumps during the last 2 decades.

Based on studies conducted by EVADIAC, the main benefits of CIPII therapy were described as follows:⁵

- Improvement of metabolic control and of glycemic stability
- Reduction of severe hypoglycaemia frequency (before CIPII: 15 per 100 patient-years, during CIPII: 2.5 per 100 patient-years).

A significant reduction of severe hypoglycaemia (less than half the frequency compared to CSII) was also found when using the Accu-Chek DiaPort system for CIPII therapy.⁶ While port systems and implantable pumps were not compared in randomised controlled trials, similar metabolic advantages may be attributed to both methods of CIPII therapy.

In a randomised cross-over study in 24 patients with type I diabetes and moderate or poor glycemic control on previous intensive subcutaneous insulin therapy, improved glycemic control was achieved with the use of CIPII⁷, thus confirming earlier observations:

- 0.8 % decrease in HbA_{1c}
- 11 % increase in the time spent in euglycemia.

Lassmann-Vague et al. attributed the observed improved glycemic stability to the intraperitoneal route of administration.⁸ Most insulin infused in the peritoneal cavity enters the portal circulation⁹, reaching first the liver, producing a porto-systemic insulin gradient closer to the physiology between the portal vein and the peripheral circulation. With a more rapid onset of action, faster reaching of peaks, faster return to baseline levels, plasma insulin kinetics are closer to physiology with CIPII than with CSII.

Treatment with an implantable pump has also some disadvantages, which can be avoided by using the Accu-Chek DiaPort system:¹⁰

- Pump explantation is required every few years for battery replacement
- Extra visits are required for re-filling the insulin reservoir of the pump every 2 to 3 months
- The large size of the implantable pump may not be suitable for children or adults with limited areas of subcutaneous tissue.

⁴ Broussolle C., et al., *French multicentre experience of implantable insulin pumps*. The Lancet 1994;343:514-515.

Renard E., Schaepelynck-Bélicar P., EVADIAC Group, *Implantable Insulin Pumps. A Position Statement About Their Clinical Use*. Diabetes Metab, 2007 Apr;33(2):158-66.

⁵ Lassmann-Vague V., Schaepelynck-Bélicar P., *The EVADIAC experience*. Diabetes Research and Clinical Practice, 2006 Dez;74:164-168.

Hanaire-Broutin H., et al., *Feasibility of Intraperitoneal Insulin Therapy With Programmable Implantable Pumps in IDDM. A Multicenter Study*. Diabetes Care, 1995 Mar;18(3):388-92.

⁶ Liebl A., et al., *A Reduction in Severe Hypoglycaemia in Type 1 Diabetes in a Randomized Crossover Study of Continuous Intraperitoneal Compared With Subcutaneous Insulin Infusion*. Diabetes Obes Metab, 2009 Nov;11(11):1001-8.

⁷ Logtenberg S.J., et al., *Improved Glycemic Control With Intraperitoneal Versus Subcutaneous Insulin in Type 1 Diabetes: A Randomized Controlled Trial*. Diabetes Care, 2009 Aug;32(8):1372-7.

⁸ Lassmann-Vague V., Schaepelynck-Bélicar P., *The EVADIAC experience*. Diabetes Research and Clinical Practice, 2006 Dez;74:164-168.

⁹ Broussolle C., et al., *French multicentre experience of implantable insulin pumps*. The Lancet 1994;343:514-515.

¹⁰ Medtronic MiniMed Inc., *Implantable Insulin Pump MiniMed(TM) 2007 Implantable Pump Patient Manual* (size of pump)

Haveman J.W., Logtenberg S.J., Kleefstra N., Groenier K.H., Bilo H.J., Blomme A.M., *Surgical Aspects and Complications of Continuous Intraperitoneal Insulin Infusion With an Implantable Pump*. Langenbecks Arch Surg, 2010 Jan;395(1):65-71.

Patient experiences with the port

Dr. Andreas Liebl, Diabetologist, Germany, over 20 years experience with CIPII, reported in 2015 on his experience with the previous generation of the port:

“With the port we have 2 important advantages: one is an excellent pharmaco kinetics and pharmaco dynamics of the insulin infusion. That is the insulin which is infused through the port acts nearly as fast as intravenously given insulin – so it’s a very quick onset and a quick offset of the insulin and I can react very quickly to blood glucose fluctuation.

And the second advantage is the fact that the distribution of the insulin in the human body is like normal physiology – which means that we have a high concentration of insulin around the liver in front of the liver and a lower concentration of insulin in the periphery.

And currently we have approximately 25 patients on port which are treated in our clinic. A patient who is on port and the port is working well, never will change his therapy. A patient who is on a port loves his port and also the doctors love a port. Because as long as the port is working properly the quality of metabolic control is better than with any other therapy we know currently. Therefore a patient never voluntarily will change his therapy.

Several things have been improved with the new port. One thing is that the number of parts which we have, has been reduced. So it’s much easier to handle for example the exchange of the membrane. It is very easy to do. Another thing which we have improved is the biocompatibility of the materials; so the materials used are very friendly to the patients’ tissue. Especially the catheter which is in the abdomen is very soft and is very tissue friendly and so we believe that we get fewer problems with adhesion in the abdomen. Also the new port is surrounded with a special tissue which allows the tissue of the body to grow into it and therefore we also believe that the number of infections will be reduced a lot with this new design. Then the surgical procedure is extremely simple and it takes just something like 15 to 20 minutes to implant the port.

So after the port has been implanted the patients are advised to rest as much as possible for 1 or 2 days. The in-growth of the material and the healing is very fast with the new port. The patient can go back to his normal activities very quickly. And there are virtually no limitations what the patient can do with the port. Many things can be done even better because the control of the blood glucose is much better with the port. So it’s more flexible and he has more options to do.”

2 System Components

The port system consists of an infusion set with ball cannula, a membrane, a fixation disc, a metal (titanium) port body, and a catheter that extends into the abdomen.

All materials used in the port system are biocompatible and hypoallergenic. They do not contain latex or PVC. Therefore, the risk of the port being rejected by the body is very remote.

The previous generation, on the market since 2000, has been replaced by a newly designed and improved port system. The improved features are:

- A special felt material is wrapped around the stem of the port body to enable the tissue to granulate closer to the stem of the port. This acts as an infection barrier.
- The new catheter material is more flexible, moves more freely and has a softer tip. This helps to reduce the risk of adhesions or overgrowing of tissue on the catheter tip.
- The risk of occlusion is expected to be significantly reduced due to an enlarged inner diameter of the catheter and the trumpet-shaped soft tip.

2.1 Port Body

Membrane (1), port body (2), and catheter (3)

The port enables intraperitoneal insulin infusion: That is, infusion of insulin into the peritoneal cavity using an Accu-Chek Spirit Combo insulin pump and a specially designed stainless steel ball cannula infusion set.

The top of the port body (with a diameter of 9 mm (0.35 inches)) sits above the surface of the skin up to 5 mm (0.19 inches). A flower-shaped plate (with a diameter of 26 mm (1.02 inches)) is placed under the skin. The fixation disc and the polyester felt band provide stability for the port system. The insulin is infused into the intraperitoneal cavity via the ball cannula infusion set.

PRECAUTION

Risk of adverse reaction

In rare situations, the implanted system may cause adverse reactions.

WARNING

Risk of scarring

In rare situations, the port body may be rejected by the human body, which requires explantation. Scars remain after explantation of the port system.

Materials of the port body:

Component	Materials	Volume	Ratio
Base	Titanium 6 AL 4V ELI + Anodisation	630.43 mm ³ (0.038 in ³)	75.20 %
Glue	Silicone adhesive MED-1511	Traces	Traces
Primer	MED6-161	Traces	Traces
Polyester felt	Polyester felt 6077 (Dacron)	207.94 mm ³ (0.013 in ³)	24.80 %

Expected lifetime: 5 years

It is recommended not to explant the port after 5 years if the port is still functional.

2.2 Membrane

The membrane seals the access to the peritoneal cavity. When the ball cannula is removed, the port is automatically sealed tight by the silicone membrane. No water or solution can go through the port.

The membrane must be changed every 6 months to prevent leakage. This is a virtually pain-free procedure, but must be performed at a center of excellence by a trained healthcare professional or surgeon under sterile conditions.

Materials of the membrane:

Component	Materials	Volume	Ratio
Lid	Titanium 6 AL 4V ELI + Anodisation	119.89 mm ³ (0.007 in ³)	62.67 %
O-ring	SILOPREN LSR 4040 (40 shore A)	9.19 mm ³ (0.00056 in ³)	4.80 %
Valve	SILOPREN LSR 4060 (60 shore A)	62.02 mm ³ (0.004 in ³)	32.50 %

Replacement interval: At least every 6 months

2.3 Catheter

The catheter enables the infusion of insulin from the external insulin pump into the peritoneal cavity. The catheter is implanted into the abdomen and is not visible outside of the body. The catheter is flexible, and has a trumpet-shaped soft tip so that it can move around in the abdominal cavity. This free movement helps to reduce the risk of adhesions or overgrowing of tissue on the catheter tip.

3 different lengths are available, 9 cm, 15 cm and 25 cm (3.5 inches, 5.9 inches and 9.8 inches). The surgeon decides which length is suitable for the patient.

The catheter must be replaced if it is blocked or overgrown. The replacement must be performed at a center of excellence by a trained surgeon under sterile conditions.


Materials of the catheter:

Component	Materials	Volume	Ratio
Catheter	Tecothane TT2074A-B20	238.05 mm ³ (0.0145 in ³)	28.30 %
Catheter	Tecothane TT1074A	521.28 mm ³ (0.0318 in ³)	61.90 %
Catheter	PTFE F polyflon	30.28 mm ³ (0.0184 in ³)	3.60 %
Hub	Tecothane TT1065D	52.71 mm ³ (0.0032 in ³)	6.30 %

Expected lifetime: 5 years

It is recommended not to explant the catheter after 5 years if the catheter is still functional.

2.4 Fixation Disc

Fixation disc (1) and disc applicator (2)

The fixation disc provides stability and flexibility. This is important because the port is not sutured into the subcutaneous tissue. Instead, the port takes hold through progressive in-growth/granulation of the subcutaneous tissue.

The implanted system with the fixation disc

Even when the healing process is complete, the fixation disc continues to stabilise the port, absorb impacts, and is also helpful in holding the port when connecting and disconnecting the infusion set.

The Accu-Chek DiaPort fixation disc is fitted to the port after implantation. This helps to reduce stress on the tissue surrounding the port site by limiting movement within the implantation site/pocket. Additionally, it helps to prevent the port from sinking or being overgrown by skin. The disc applicator is used for fitting and removing the fixation disc.

Patients must wear the fixation disc at all times. Change the fixation disc when changing the cartridge of the insulin pump or at least every 6 days. For more information, see all warnings and cautions, hygiene requirements, product use or care and maintenance recommendations for the Accu-Chek DiaPort system in the *Handbook for Port Users* and in the instructions for use that come with the product.

Materials of the fixation disc and disc applicator:

Component	Materials	Volume	Ratio
Fixation disc	PP BORMED RF825 MO	n/a	n/a
Disc applicator	Compound based on PP 1013 H1	n/a	n/a
Disc applicator	Blue colorant Lifocolor PP (0.5 %)	n/a	n/a

Replacement interval: At least every 6 days

2.5 Infusion Set

Infusion set

The infusion set connects the system with an external Accu-Chek Spirit Combo insulin pump. When patients use the infusion set for the first time, a healthcare professional must be present.

The infusion set has been specially designed for use with the port, and has a ball cannula that is inserted into the port body membrane, which enables the connection with the catheter. This solution avoids needle-stick injuries, and is virtually pain free. After disconnection of the infusion set, the membrane automatically seals the access to the peritoneal cavity.


The infusion set with the handling aid attached to the insulin pump

Replacement interval: At least every 6 days

3 Implantation

3.1 Implantation Set

The implantation set is composed of 2 sets.

The Accu-Chek DiaPort implantation set: Part I
Disc applicator (1), screw-in screwdriver with a pre-mounted membrane (2), port body with pre-mounted catheter (15 cm (5.9 inches)) (3), fork wrench (4), 2 catheters (9 cm (3.5 inches) and 25 cm (9.8 inches)) (5) and fixation disc (6)

The Accu-Chek DiaPort implantation set: Part II
Layer scalpel (1), syringe (2), infusion set (3), skin punch (4) and punch plate (5)

3.2 Surgical Procedure

The port is implanted under general anesthesia during a minor operation that takes approximately 20 to 30 minutes. Discuss the placement of the port, the implantation, and the anesthesia in advance with the diabetologist, surgeon, and anesthetist.

The success of the implantation depends on avoiding an infection. Before surgery, it is recommended to carefully examine the abdomen of the patient. There should be no sign of inflamed or infected tissue, as this could lead to complications after implantation. If necessary, the implantation should be postponed until the area has healed.

The operation takes place in a hospital that is associated with the center of excellence. For people with type I diabetes, strict glycemic control is recommended during surgery. It is recommended that patients remain in hospital for at least one-and-a-half days. There may be slight to moderate pain following the operation, but this can be managed with analgesics.

WARNING

Risk of abdominal pain

The catheter of the Accu-Chek DiaPort system is placed in the abdominal cavity. Although the catheter is made of a flexible material, it can cause abdominal pain.

The implant card and the implantation set sticker are dispensed with the Accu-Chek DiaPort implantation set. The filled in implant card has to be given to the port user before leaving the hospital.

This set contains pointed or sharp-edged components. Dispose of pointed or sharp-edged components according to local regulations, in such a way that no third party can be injured.

For detailed handling instructions, refer to the instructions for use provided with the product.

3.3 Implantation Site

The port can be implanted wherever the patient prefers in the lower abdominal region as shown by the green shaded area in the pictures, with a few exceptions. Do not position the port on or near pressure points.

For example, the port must **not** be positioned in:




- the region close to the navel.
- the region crossed by a car safety belt.
- the region under or near the trouser belt.
- the region around the rib-arch.
- a place where it can bump against the edge of a table (when seated).

Right-handed patients may prefer to have the port on the left side, and vice versa.

It is important that the patient's daily activities are not negatively affected by the positioning of the port. Therefore, carefully discuss with him or her the different activities to be undertaken in school, work or leisure time, especially sporting activities.

Make sure to discuss a 2nd surgical placement option in case there is an issue during surgery, preventing the port from being placed at the original preferred site.

Select the implantation site in such a way that obese patients can see the port during handling.


3.4 What Will the Implanted Port Look Like?

The operation leaves a scar which is about 4 cm (1.57 inches) long. The exact size depends on the surgical requirements during the operation (for example, unexpected complications).

The top of the port (diameter of 9 mm (0.35 inches)) protrudes above the surface of the skin by up to 5 mm (0.19 inches).

The picture shows the port with the fixation disc and the infusion set several weeks after surgery.

The implanted Accu-Chek DiaPort system

4 Postoperative Care

Normally, the port can be used immediately following implantation, but this decision is under the responsibility of the diabetologist/endocrinologist.

After the implantation, the healing process includes the healing of the surgical wound and of the implantation pocket in which the port is in-growing. The scar left by the operation, needs about 10 days to heal. Healing of the tissue surrounding port implantation site normally takes between 4 to 6 weeks, if no infection or other problems occur.

Adequate rest

During the first few days after implantation, it is recommended that the port user rests as much as possible and remains in hospital for at least one-and-a-half days. This supports the healing of the wound and enables you to make any necessary adjustments to your patients diabetes management plan as the patient's body starts to adjust to CIPII therapy.

It is very important that your patient does not undertake activities that affect the healing of the wound. It is strongly recommended not to lift or carry heavy objects during the healing time.

Dressing the port site and the surgical wound

The success of the port therapy depends on how well the tissue grows around the port, and how clean the port is kept. Therefore, you must instruct your patient on how to keep the port body area as clean as possible.

The port and the surgical wound should be cleaned at least once a day and dressed with a dry, sterile dressing (padded if necessary) until the wound has healed completely. It is recommended to only use sterile/clean water, or a physiological solution for cleaning the port body. The cleaning should be done without any pressure or rubbing using normal aseptic techniques.

After showering, and before putting a new dressing on, the skin should be dried carefully using a towel.

It must be emphasised to your patient that they must not graze or pick the scab growing on the wound around the port. This is part of the healing process.


Other Precautions

For more information on precautions to be taken by patients, see the *Handbook for Port Users*.

5 Therapy

5.1 Insulin Kinetics With CIPII Therapy

A faster onset of insulin action, and a shorter duration of action were observed after intraperitoneal insulin infusion.¹¹

Insulin concentrations during intraperitoneal insulin infusion and subcutaneous insulin injection.¹¹

While the insulin kinetics during CIPII therapy more closely mimic the physiologic pattern of insulin secretion, and contribute to improved glycemic control and a lower frequency of severe hypoglycaemia¹², the differences between CIPII and CSII need to be taken into account during insulin pump therapy.

WARNING

Risk of hyperglycaemia or hypoglycaemia

When starting CIPII therapy, the patient's insulin needs are likely to change. Speak to the patient about the adjustments to his or her therapy parameters.

The patient's therapy settings, especially the bolus advice settings, need to be carefully looked at and if necessary adapted when initiating CIPII therapy to take into account differences in insulin kinetics and glucose profiles, with regard to the faster onset and shorter duration of insulin action.

Summary

The kinetics of the insulin infused using the intraperitoneal pathway is much faster than with subcutaneous insulin infusion. This means that some therapy parameters like meal factors, and more especially for bolus calculation, will need to be adapted. The use of the different bolus patterns needs to be adapted to the much faster kinetics of insulin.

¹¹ Schade D.S., et al., *Normalization of plasma insulin profiles with intraperitoneal insulin infusion in diabetic man*. Diabetologia 1980; 19(1):35-39.

¹² Nathan D.M., et al., *Postprandial Insulin Profiles With Implantable Pump Therapy May Explain Decreased Frequency of Severe Hypoglycemia, Compared With Intensive Subcutaneous Regimens, in Insulin-Dependent Diabetes Mellitus Patients*. Am J Med, 1996 Apr;100(4):412-7.

5.2 Switching from CSII to CIPII

The switch from subcutaneous insulin pump therapy (CSII) to intraperitoneal insulin therapy (CIPII) can take place immediately after the implantation of the port. The insulin pump can be connected to the port as soon as the implantation has been completed.

Clinical experience has shown that rinsing the catheter after implantation for several hours with physiological saline solution (1 mL per hour, for up to 24 hours) helps to avoid infections, blockages in the catheter, and insulin induced irritation of the peritoneum.

Depending on the indications for using the port, you may need to adapt the insulin requirement:

- In case of confirmed subcutaneous insulin resistance, the total daily dose will usually need to be significantly reduced as the insulin may act much more effectively. Closely monitor the glucose values during the first days until the dose has been adapted to the new needs of the patient.
- For the other indications, a general reduction of 10 % – 15 % of the total daily dose is recommended at the beginning of the therapy and has to be adapted on an ongoing basis according to the needs of the patient. In general the insulin doses required for intraperitoneal therapy are similar to those required in intravenous therapy.

The following points for the use of bolus calculation have to be considered:

- The parameters for bolus calculation need to be re-defined. As the duration of insulin action is shorter under CIPII therapy, specific parameters have to be reviewed and adapted.
- Acting time and offset time usually need to be reduced. Please refer to the corresponding instructions for use to carry out these adjustments.

Which insulin can be used?

Not all types of insulin are suitable for CIPII therapy. **Use only U100 short-acting regular human insulin.** We recommend using **Insuman® Infusat** because it showed a low tendency towards occlusions after several weeks.¹³

Preparation and connection of the Accu-Chek Spirit Combo insulin pump with the Accu-Chek DiaPort system

Using an Accu-Chek insulin pump with the port is very similar to using a pump for subcutaneous insulin infusion. The insulin pump can be used in the same way.

WARNING

Risk of hyperglycaemia

The use of a non-compatible insulin pump can lead to insufficient insulin delivery. Use the Accu-Chek DiaPort system only with an Accu-Chek Spirit Combo insulin pump.

Risk of hypoglycaemia

When filling the infusion set with insulin, the insulin pump delivers a considerable insulin amount at once. Do not fill the infusion set while it is connected to the body. You may deliver an uncontrolled insulin amount.

¹³ Liebl A., et al., *Evaluation of the New ACCU-CHEK® DIAPORT, a Port System for Continuous Intraperitoneal Insulin Infusion, in Patients With Type 1 Diabetes: First 6-Month Results.* Diabetes, 2013; 62: A247-A248.

Before starting therapy with the port, change the insulin cartridge of the Accu-Chek Spirit Combo insulin pump. Make sure the infusion set is disconnected from your patient. Fill the infusion set according to the instructions in the User's Manual of the insulin pump. The fill amount of the infusion set depends on its length:

- 70 cm (27.56 inches) infusion set: 13 U (U100)
- 100 cm (39.37 inches) infusion set: 17 U (U100)

When connecting the pump for the first time after the implantation, the implanted catheter must also be filled. The fill amount of the catheter depends on the length of the implanted catheter:

- 9 cm (3.5 inches) catheter: 3.5 U (U100)
- 15 cm (5.9 inches) catheter: 6 U (U100)
- 25 cm (9.8 inches) catheter: 10 U (U100)

When delivering the first insulin dose through CIPII, a healthcare professional must be present. It is strongly recommended to closely observe the insulin regime using more frequent blood glucose monitoring.

How long can the insulin pump remain disconnected?

Compared to CSII therapy, no subcutaneous insulin depot is present. This means that a few minutes after the insulin pump has been disconnected, very little insulin is still available in the body.

During the disconnection time, the same rules as for CSII therapy can be followed with respect to the insulin kinetics.

If the patient disconnects the pump for a short time, for example, when taking a shower, it is recommended **not** to put the pump in STOP mode. The insulin dripping from the infusion set keeps the tip clean and reduces the risk of air bubbles in the tubing. Patients can use the handling aid to attach the infusion set on the luer-lock connector of the pump. Do not touch the ball cannula directly after removal of the protective cap.

5.3 Hypoglycaemia and Hyperglycaemia

Mild hypoglycaemia

Mild hypoglycaemia can be treated the same way as with CSII therapy.

Severe hypoglycaemia

In the event of severe hypoglycemia, interrupt the insulin delivery immediately. This is done by either stopping the insulin pump or by removing the ball cannula from the port. Additionally, patients should eat carbohydrates as recommended by their healthcare professional and control their blood glucose regularly.

Hyperglycaemia and/or diabetic ketoacidosis

If the patient experiences symptoms of hyperglycaemia and/or diabetic ketoacidosis, he or she should follow the treatment plan as recommended by their healthcare professional.

First, eliminate all reasons that may be related to the insulin pump or the infusion set. Make sure the infusion set is properly connected and that there are no air bubbles in the tubing. If you think that the infusion set may not be working properly, change it immediately.

If changing the infusion set does not help, check if the catheter is blocked. **Carefully** flush the catheter with physiological saline solution using the 3 mL syringe. During this procedure, the content of the catheter will be infused.

WARNING

Risk of hypoglycaemia

When flushing the catheter or when inserting the guide wire during catheter replacement, the content of the catheter is delivered to the patient's body. The amount of insulin depends on the length of the implanted catheter:

- 9 cm (3.5 inches) catheter: 3.5 U (U100)
- 15 cm (5.9 inches) catheter: 6 U (U100)
- 25 cm (9.8 inches) catheter: 10 U (U100)


Closely monitor your patient's blood glucose during and after this procedure.

Control your patient's blood glucose during this procedure, as the risk of hypoglycaemia is significantly increased. Continue the tight monitoring for at least one week to check if there is any improvement.

In case of no improvement, check whether the catheter is overgrown. To check whether the catheter is overgrown, a small amount of contrast media can be used to see if there is any insulin flooding back outside of the catheter. Backflow of insulin can be a sign of adhesions or encapsulation, which can be carefully removed via laparoscopy. If the catheter is occluded, a catheter replacement is necessary.

5.4 Catheter Replacement

The catheter replacement set contains the following components:

The Accu-Chek DiaPort catheter replacement set

Guide-wire (1), screw-in screwdriver with pre-mounted membrane (2), catheter extractor (3), screw-off screwdriver (4), fork wrench (5), catheters in three different lengths (9 cm, 15 cm and 25 cm (3.5 inches, 5.9 inches and 9.8 inches)) (6), fixation disc (7) and disc applicator (8)

The catheter is connected to the infusion set via the membrane and transfers the insulin from the external pump into the peritoneal cavity. 3 catheter lengths are available:

- 9 cm (3.5 inches) for children or people with a low body mass index (BMI)
- 15 cm (5.9 inches) for people with a medium BMI
- 25 cm (9.8 inches) for people with a high BMI

The Accu-Chek DiaPort catheter replacement set is a surgical toolkit for changing the catheter of the Accu-Chek DiaPort system. The product is intended only for use by healthcare professionals trained in using the Accu-Chek DiaPort system, using the tools provided with the set. Catheter replacement must be performed in a clean environment under sterile conditions.

WARNING

Risk of injury

If you modify the system or if you use tools that were not provided with the set, you may injure internal organs.

Do not modify the system. Use only the tools provided with the set and adhere to the instructions for use.

The catheter should only be replaced if it is blocked or overgrown by tissue in the peritoneal cavity. Salvage of a blocked catheter can either be achieved by a laparoscopic approach under anesthesia to free the peritoneal tip of the catheter, or by an out-patient procedure without anesthesia.

 **PRECAUTION**

Risk of interruption of therapy

If the catheter is blocked or overgrown, a laparoscopic approach may be required to free the tip of the catheter.

Make sure that the patient has an alternative insulin therapy available.

This set contains pointed or sharp-edged components. Dispose of pointed or sharp-edged components according to local regulations, in such a way that no third party can be injured.

 **WARNING**

Risk of hypoglycaemia

When flushing the catheter or when inserting the guide wire during catheter replacement, the content of the catheter is delivered to the patient's body. The amount of insulin depends on the length of the implanted catheter:


- 9 cm (3.5 inches) catheter: 3.5 U (U100)
- 15 cm (5.9 inches) catheter: 6 U (U100)
- 25 cm (9.8 inches) catheter: 10 U (U100)

Closely monitor your patient's blood glucose during and after this procedure.

For detailed handling instructions, refer to the instructions for use provided with the product.

5.5 Membrane Replacement

The membrane replacement set contains the following components:

*The Accu-Chek DiaPort membrane replacement set
Disc applicator (1), screw-in screwdriver with pre-mounted membrane (2),
screw-off screwdriver (3), fork wrench (4), and fixation disc (5)*

The membrane is a component of the port body that needs to be exchanged every 6 months to prevent leakages. The ball cannula tip connects through the membrane into the hub of the catheter. After disconnection of the infusion set, the membrane automatically seals the access to the peritoneal cavity.

The Accu-Chek DiaPort membrane replacement set enables the change of the Accu-Chek DiaPort membrane. The product is intended only for use by healthcare professionals trained in the Accu-Chek DiaPort system, using the tools provided with the set. Membrane replacement must be performed in a clean environment under sterile conditions.

Prior to membrane/catheter replacement, it may be advisable to take a 30 minutes warm bath to dissolve leftovers or crystallisation from insulin.

WARNING

Risk of hyperglycaemia

A dislodged catheter hub can lead to insufficient insulin delivery.

Make sure that the catheter hub is correctly aligned during and after the change of the membrane. The notch of the catheter hub fits into the port body in a defined direction only.

This set contains pointed or sharp-edged components. Dispose of pointed or sharp-edged components according to local regulations, in such a way that no third party can be injured.

For detailed handling instructions, refer to the instructions for use provided with the product.

6 Daily Care

Prevention of infection and inflammation

It is important to avoid any infection or inflammation around the port. This is one of the most frequent causes for explantation of the port.

To prevent infection and inflammation of the implantation site, urge the patient to follow these guidelines:

- Always handle the port with clean dry hands.
- Check the site of the port once a day, to make sure that everything is alright.
- Do not place cream around the port area.
- Clean the port and the surrounding area with sterile or clean water, or with a physiological saline solution. In general it is not necessary to use a skin disinfectant. Do not use solutions containing alcohol, as these can be too aggressive for the skin.
- Carefully dry the skin around the port with a clean towel, without pressure or rubbing.
- Always cover the area with an adhesive dressing if there are risks of soiling or infection. Make sure the adhesive dressing is large enough so that it does not stick to the membrane or the port body.

In the event of inflammation or infection

As soon as any swelling, pain, oozing, or reddening is visible at the port site, experienced healthcare professionals from the certified Accu-Chek DiaPort centers of excellence recommend to use a skin disinfectant (for example, containing octenidine dehydrochloride) in addition to the normal cleaning procedure. It may be advisable to start antibiotic therapy in addition. Start a culture to identify the pathogenic agent. If necessary, start intravenous antibiotics therapy.

If the treatment of the infection is not successful, the port must be removed and the treatment with antibiotics continued until the infection has disappeared.

WARNING

Risk of scarring

In some situations, the port body must be explanted, for example, if medical complications occur. Scars remain after explantation of the port system.

However, once the infection is cleared, a new port can be implanted at a new site.

7 Medical Conditions not Directly Related to the Port

Pregnancy

Pregnancies with the port are known without any complications. No removal of the port was necessary. The therapy type was not changed during the pregnancy, the HbA_{1c} was 5.6 – 5.8 %, and no severe hypoglycaemia was recorded.¹⁴

If one of your patients plans to become pregnant, it is recommended to increase the frequency of their visits and to monitor the skin around the port more often. Apart from this, no additional steps are required to be taken, beyond those for continuous subcutaneous insulin infusion therapy.

Medical Imagery

The catheter is designed in a way that it appears on X-rays and Magnetic Resonance Imaging pictures.

However, the insulin pump should be removed from the room before X-ray or Magnetic Resonance Imaging examination to avoid any electronic errors. For more information, please refer to the warnings and precautions of the insulin pump manual.

Other imaging methods like nuclear medicine or ultrasound scans should not affect, nor be affected by, the port.

Consultations with your patient

As with continuous subcutaneous insulin therapy, the patient is encouraged to visit you at least every 3 months, or as necessary. Note that membrane replacement is required every 6 months. If the port is working correctly, no additional visits are required. Your patient is also encouraged to inform you whenever questions arise, or if complications or difficulties occur.

Hospitalisation

If the patient has to be hospitalised, even for non-diabetes therapy-related reasons, he/she is urged to discuss this with you.

Emergency situations

If the patient is admitted to hospital in an emergency (not related to the port) or as a result of an accident, it is important that the hospital staff is aware of the port system including the intraperitoneal catheter, in case surgery is required. It is recommended that patients always carry the Accu-Chek DiaPort implant card received after the implantation. It contains information on the implantation, the surgeon and the institution that performed the implantation, so that the medical staff can get additional information if required.

¹⁴ Fonseca V.A., et al., *Diabetic Pregnancy Managed With Intraperitoneal Insulin*. Diabet Med, Jan-Feb 1987;4(1):74-6.

Fonseca V.A., et al., *Successful Pregnancy in Diabetic Controlled With Intraperitoneal Insulin*. Diabetes Care, Jul-Aug 1987;10(4):541-2.

Schnell O., et al., *A case of diabetic pregnancy controlled with a percutaneous access device for intraperitoneal insulin infusion*. Diabetes Care. 1994;17(11):1354-1355

8 Daily Life

Day-to-day activities and therapy with an Accu-Chek Spirit Combo insulin pump should not be affected by the port system. Below is a summary of the advices contained in the *Handbook for Port Users*.

Continuing work

Except for manual labor or work which may involve risks associated with contact or pressure on the port system, the port should have little impact on your patients professional activities.

It is particularly important, in the weeks following the operation, that patients avoid contact or pressure on the port, so that the wound can heal properly.

Driving a vehicle

When choosing the site for the implantation, take into account that patients may have to drive a vehicle. Therefore, it is crucial to position the port so that the safety belt does not interfere with the port or the wearing of the insulin pump.

Physical work

Once the healing process of the subcutaneous pocket and tissue surrounding the port body is complete, generally after 4 to 6 weeks, there should be no problem with lifting or carrying weights.

Your patient should consult you prior to lifting heavy weights or undertaking physical work for the first time post implantation.

Sport

Participating in sporting activities requires some care and common sense. The first rule is that the healing process must be complete, with no significant inflammation. In addition, the port must not be subjected to strong pressure or strain.

Activities such as hiking or cycling generally do not present difficulties. However, it is not recommended to participate in contact sports with risk of injury, such as football or boxing.

Swimming

Swimming is permitted, provided that patients only swim in pools with chlorinated water, clean lake water, or in the sea. This is in order to help prevent infections.

It is recommended to cover the port with a waterproof adhesive dressing before swimming. The adhesive dressing must be large enough so that it does not stick to the membrane or the port body. After swimming, the waterproof adhesive dressing should be removed and the port should be carefully dried with a clean towel.

Sunbathing/Solarium

The port should not be exposed to direct sunlight for any length of time, as the port body can become very hot and burn the skin around the port. This may affect the stability of the port and can be a cause of infection.

WARNING

Risk of skin burns

If the port body is exposed to high temperatures, for example, in a sauna, skin burns can occur. Do not expose the port body to high temperatures.

However, if the patient wishes to sunbath, he or she should be carefully trained to cover the port, and a large area around it, with an adhesive dressing to protect the skin. The dressing must be placed in such a way that the port does not heat up to a point where it could burn the surrounding skin.

Showering

During the 3 days after implantation, it is not recommended to take a shower. Afterwards, when showering, it is recommended that patients take this opportunity to clean the port and the surrounding skin at least once a day.

During the healing period: After the shower, the port and the surgical wound must be dressed with a dry, sterile dressing (padded if necessary) until the wound has healed completely.

After the healing period: After the shower, patients should always dry the port very carefully with a clean towel, without pressure or rubbing.

Sexual intercourse

There are no additional requirements.

Sauna

Sauna is not recommended, as the port body can become very hot and burn the surrounding skin. Sauna also increases the risk of infection.

Sleeping on the abdomen

Sleeping on the abdomen seems not to be an issue as reported by patients.

Flying

Travel by aircraft does not present problems that are specific to the port. The port is detected by security portals because it is made with titanium, which is already used for many implanted devices.

With the help of the implant card, patients can explain what the port is. Normal precautions for the insulin pump still apply. For more information, refer to the User's Manual of the insulin pump.

Travelling

If patients have been using continuous subcutaneous insulin infusion (CSII) until now, they can manage travelling in a similar way.

However, situations may arise in which patients may need specialised medical attention: For example, if the catheter is blocked. Therefore, patients are advised to carry the implant card at all times. Furthermore, patients should have a backup plan and should take extra supplies with them, including insulin pens or needles to switch from port therapy in case of a malfunction of the port.

Glossary

Allergy	Over-sensitive reaction to certain materials
Catheter	A flexible tube made of synthetic material that delivers the insulin into the abdominal (peritoneal) cavity.
CIPII	Continuous intraperitoneal insulin infusion
CSII	Continuous subcutaneous insulin infusion
Disc applicator	A tool which is supporting the connection and disconnection of the fixation disc to the body of the port.
Fixation disc	A plastic ring that stabilises the port. It should be kept in place at all times.
Handling aid	An aid that helps to connect or to disconnect the infusion set from the port.
Hyperglycaemia	High blood glucose level
Hypoglycaemia	Low blood glucose level
Infusion set with ball cannula	A special connecting circuit between the insulin pump and the port
Insulin resistance	A sharp drop in the therapeutic effect of insulin, or its failure to work due to insensitivity of the body to insulin
Lipodystrophy	A pathological increase or decrease in the fatty tissue
Percutaneous port system	A port system, where some parts are inside and outside of the body.
Peritoneal	Concerning the space in the abdomen lined with the peritoneum
Physiological effect	An effect that corresponds to the normal functioning of an organism.

Safety loop

The tubing of the infusion set is taped to the skin near the infusion site in the shape of a loop to avoid dislodgement of the cannula due to accidental dragging on the tubing.

Subcutaneous tissue

Tissue under the surface of the skin

Symbol Descriptions

	Consult instructions for use or consult electronic instructions for use
	Caution, refer to safety-related notes in the instructions for use accompanying this product.
	Temperature limit
	Keep away from sunlight
	Keep dry
	Use by (YYYY-MM-DD)
	Do not use if package is damaged and consult instructions for use
	Sterilised using ethylene oxide
	Single sterile barrier system with protective packaging outside
	Double sterile barrier system with protective packaging outside
	Use only once
	Date of manufacture
	Medical device
	Manufacturer
	Unique device identifier
	Catalogue number
	Batch code
	Complies with the provisions of the applicable EU Legislation

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Tel: 0800-93626 (Accu-Chek Service)
www.accu-chek.be

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Komponenten des Accu-Chek DiaPort Systems:

<p>1</p>  A blue, dome-shaped plastic component with horizontal ridges and a small protrusion at the bottom. The word "ACCUCHEK" is embossed on the side.	<p>Platzierungshilfe</p>
<p>2</p>  A clear plastic infusion set with a curved tube and a small metal ball cannula at the end.	<p>Infusionsset mit Kugelkanüle</p>
<p>3</p>  A blue, flexible plastic applicator with a curved, ribbed design.	<p>Scheibenappikator</p>
<p>4</p>  A clear plastic circular disc with a central hole and a small protrusion.	<p>Fixationsscheibe</p>

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

Riveline J.P., et al., *Subcutaneous Insulin Resistance Successfully Circumvented on Long Term by Peritoneal Insulin Delivery From an Implantable Pump in Four Diabetic Patients*. Diabetes Metab, 2005 Nov;31(5):496-8.

Baillot-Rudoni S., et al., *Implantable Pump Therapy Restores Metabolic Control and Quality of Life in Type 1 Diabetic Patients With Buschke's Nonsystemic Scleroderma*. Diabetes Care, 2006 Jul;29(7):1710.

Hilgard D., Laengler A., *Severe subcutaneous insulin infusion resistance in a 13 year old child - a case-report of successfully performed therapy with a diaport™-system*. Pediatric Diabetes, 2009; 10(Suppl 11):37.

Noud M.N., et al., *Benefits of intra-peritoneal insulin administration in a child with severe insulin-induced lipoatrophy*. Pediatric Diabetes, 2009;10 (Suppl 11):100.

Lee S.I., Narendran P., *Intraperitoneal Insulin Therapy for a Patient With Type 1 Diabetes With Insulin Injection Site Inflammation*. BMJ Case Reports, 2014 Aug 21.

Garcia-Verdugo R., Erbach M., Schnell O., *A New Optimized Percutaneous Access System for CIPII*. J Diabetes Sci Technol. 2017 Jul;11(4):814-821..

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Componenti del sistema Accu-Chek DiaPort:

<p>1</p>  A blue, dome-shaped plastic guide device with a small protrusion at the bottom. The word "ACCUCHEK" is embossed on the front.	<p>Dispositivo di guida</p>
<p>2</p>  A clear plastic infusion set with a thin tube and a small metal cannula at the end.	<p>Set d'infusione con cannula a sfera</p>
<p>3</p>  A blue, flexible plastic applicator with a curved shape and a small notch.	<p>Applicatore del disco</p>
<p>4</p>  A clear plastic circular disc with a central hole and a small protrusion.	<p>Disco di fissaggio</p>

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
¹ Schade D.S., Duckworth W.C., *In Search of the Subcutaneous-Insulin-Resistance Syndrome*. N Engl J Med, 1986 Jul 17;315(3):147-53.
Soudan B., et al., *Extreme Subcutaneous Insulin Resistance: A Misunderstood Syndrome*. Diabetes Metab, 2003 Nov;29(5):539-46.
Riveline J.P., et al., *Subcutaneous Insulin Resistance Successfully Circumvented on Long Term by Peritoneal Insulin Delivery From an Implantable Pump in Four Diabetic Patients*. Diabetes Metab, 2005 Nov;31(5):496-8.
Baillot-Rudoni S., et al., *Implantable Pump Therapy Restores Metabolic Control and Quality of Life in Type 1 Diabetic Patients With Buschke's Nonsystemic Scleroderma*. Diabetes Care, 2006 Jul;29(7):1710.
Hilgard D., Laengler A., *Severe subcutaneous insulin infusion resistance in a 13 year old child - a case-report of successfully performed therapy with a diaporTM-system*. Pediatric Diabetes, 2009; 10(Suppl 11):37.
Noud M.N., et al., *Benefits of intra-peritoneal insulin administration in a child with severe insulin-induced lipoatrophy*. Pediatric Diabetes, 2009;10 (Suppl 11):100.
Lee S.I., Narendran P., *Intraperitoneal Insulin Therapy for a Patient With Type 1 Diabetes With Insulin Injection Site Inflammation*. BMJ Case Reports, 2014 Aug 21.
Garcia-Verdugo R., Erbach M., Schnell O., *A New Optimized Percutaneous Access System for CIPII*. J Diabetes Sci Technol. 2017 Jul;11(4):814-821..

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¹ Schade D.S., Duckworth W.C., *In Search of the Subcutaneous-Insulin-Resistance Syndrome*. N Engl J Med, 1986 Jul 17;315(3):147-53.
Soudan B., et al., *Extreme Subcutaneous Insulin Resistance: A Misunderstood Syndrome*. Diabetes Metab, 2003 Nov;29(5):539-46.
Riveline J.P., et al., *Subcutaneous Insulin Resistance Successfully Circumvented on Long Term by Peritoneal Insulin Delivery From an Implantable Pump in Four Diabetic Patients*. Diabetes Metab, 2005 Nov;31(5):496-8.
Baillot-Rudoni S., et al., *Implantable Pump Therapy Restores Metabolic Control and Quality of Life in Type 1 Diabetic Patients With Buschke's Nonsystemic Scleroderma*. Diabetes Care, 2006 Jul;29(7):1710.
Hilgard D., Laengler A., *Severe subcutaneous insulin infusion resistance in a 13 year old child - a case-report of successfully performed therapy with a diaport™-system*. Pediatric Diabetes, 2009; 10(Suppl 11):37.
Noud M.N., et al., *Benefits of intra-peritoneal insulin administration in a child with severe insulin-induced lipoatrophy*. Pediatric Diabetes, 2009;10 (Suppl 11):100.
Lee S.I., Narendran P., *Intraperitoneal Insulin Therapy for a Patient With Type 1 Diabetes With Insulin Injection Site Inflammation*. BMJ Case Reports, 2014 Aug 21.
Garcia-Verdugo R., Erbach M., Schnell O., *A New Optimized Percutaneous Access System for CIPII*. J Diabetes Sci Technol. 2017 Jul;11(4):814-821.

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
¹ Schade D.S., Duckworth W.C., *In Search of the Subcutaneous-Insulin-Resistance Syndrome*. N Engl J Med, 1986 Jul 17;315(3):147-53.
Soudan B., et al., *Extreme Subcutaneous Insulin Resistance: A Misunderstood Syndrome*. Diabetes Metab, 2003 Nov;29(5):539-46.
Riveline J.P., et al., *Subcutaneous Insulin Resistance Successfully Circumvented on Long Term by Peritoneal Insulin Delivery From an Implantable Pump in Four Diabetic Patients*. Diabetes Metab, 2005 Nov;31(5):496-8.
Baillot-Rudoni S., et al., *Implantable Pump Therapy Restores Metabolic Control and Quality of Life in Type 1 Diabetic Patients With Buschke's Nonsystemic Scleroderma*. Diabetes Care, 2006 Jul;29(7):1710.
Hilgard D., Laengler A., *Severe subcutaneous insulin infusion resistance in a 13 year old child - a case-report of successfully performed therapy with a diaport™-system*. Pediatric Diabetes, 2009; 10(Suppl 11):37.
Noud M.N., et al., *Benefits of intra-peritoneal insulin administration in a child with severe insulin-induced lipodystrophy*. Pediatric Diabetes, 2009;10 (Suppl 11):100.
Lee S.I., Narendran P., *Intraperitoneal Insulin Therapy for a Patient With Type 1 Diabetes With Insulin Injection Site Inflammation*. BMJ Case Reports, 2014 Aug 21.
Garcia-Verdugo R., Erbach M., Schnell O., *A New Optimized Percutaneous Access System for CIPII*. J Diabetes Sci Technol. 2017 Jul;11(4):814-821...

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© 2021 Roche Diabetes Care
ACCU-CHEK, ACCU-CHEK SPIRIT COMBO and DIAPORT are trademarks of Roche.
All other product names and trademarks are the property of their respective owners.
Roche Diabetes Care GmbH