
# Catheter replacement set

## Instructions for use

EN DE FR IT NL

**ⓘ** Consult instructions for use or consult electronic instructions for use / Gebrauchsanweisung beachten oder elektronische Gebrauchsanweisung beachten / Consulter les instructions d'utilisation ou consulter les instructions d'utilisation en format électronique / Consultare le istruzioni per l'uso in formato cartaceo o elettronico / Raadpleeg de gebruiksaanwijzing of de elektronische gebruiksaanwijzing

**⚠** Caution, refer to safety-related notes in the instructions for use accompanying this product. / Achtung, Sicherheitshinweise in der Gebrauchsanweisung dieses Produkts beachten. / Attention, se conformer aux consignes de sécurité figurant dans les instructions d'utilisation du produit. / Attenzione, osservare le avvertenze di sicurezza contenute nelle istruzioni per l'uso del prodotto. / Waarschuwing, raadpleeg de veiligheidsaanwijzingen in de gebruiksaanwijzing van het product.

Temperature limit / Temperaturbegrenzung / Limite de température / Limiti di temperatura / Temperaturbeperking

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







TI003-0 Last update: 2020-08 0044011082/6

### EN) Accu-Chek DiaPort catheter replacement set

#### Intended Use

The Accu-Chek DiaPort system is intended for continuous intraperitoneal insulin infusion (CIPI). For this purpose, it shall be used together with an external Accu-Chek Spirit Combo insulin pump. The Accu-Chek DiaPort system delivers insulin from the Accu-Chek insulin pump through the infusion set and the catheter into the peritoneal cavity.

#### WARNING Risk of hyperglycaemia

The use of a non-compatible insulin pump can lead to insufficient insulin delivery. Use the Accu-Chek DiaPort system only with an Accu-Chek Spirit Combo insulin pump.

Intended users of the Accu-Chek DiaPort system are insulin-dependent people with type 1 or type II diabetes who meet the indication for CIPI therapy. The Accu-Chek DiaPort system is only for use by patients and healthcare professionals specifically trained in the use of the individual components relative to their areas of expertise.

#### WARNING Risk of inappropriate therapy

The Accu-Chek DiaPort system is designed for intended users only. Do not use the system with patients that do not meet the indications, or where one or more contraindications apply.

The Accu-Chek DiaPort catheter replacement set is a surgical toolkit for changing the catheter of the Accu-Chek DiaPort system. The product is intended only for use by healthcare professionals trained in using the Accu-Chek DiaPort system, using the tools provided with the set. Catheter replacement must be performed in a clean environment under sterile conditions.

#### WARNING Risk of different types of harm

Modification of system components or disregard of the instructions for use may prevent the system from working as intended. This can cause different types of harm, for example, infection, hyperglycaemia, damage to the skin, or scarring due to explanation of the system. Do not modify the system and adhere to the instructions for use.

A catheter replacement is required only if the catheter is blocked, or has become overgrown by tissue in the peritoneal cavity.

Salvage of a blocked catheter can either be achieved by a laparoscopic approach under anaesthesia to free the peritoneal tip of the catheter, or by an out-patient/ambulant procedure without anaesthesia.

#### PRECAUTION Risk of interruption of therapy

If the catheter is blocked or overgrown, a laparoscopic approach may be required to free the tip of the catheter.

Make sure that the patient has an alternative insulin therapy available.

#### I Laparoscopic approach: A laparoscopic approach is required if the exchange of the catheter is due to encapsulation. The surgeon must be trained on how to perform a laparoscopic procedure.

#### II Out-patient/ambulant procedure: The procedure is described below.

#### Indications and Contraindications

See the instructions for use of the Accu-Chek DiaPort implantation set.

#### Component Overview

This is a sterile, single use product that has been sterilised with ethylene oxide.

- A Fork wrench
- B Screw-off screwdriver
- C Screw-in screwdriver with pre-mounted membrane
- D Catheter extractor
- E Catheter in 3 lengths  
Length: 9 cm 15 cm 25 cm  
Fill volume: 3.5 inches 5.9 inches 9.8 inches  
3.5 U\* 6 U\* 10 U\*
- E1 Catheter hub
- E2 Catheter tip
- F Fixation disc
- G Disc applicator
- H Guide wire
- I Sticker \*U100 insulin

#### WARNING Risk of infection

A damaged packaging may no longer be sterile or may contain a damaged component.

Before using a new component, make sure its packaging is not damaged. If the packaging is damaged, do not use its contents.

#### Risk of injury

If you modify the system or if you use tools that were not provided with the set, you may injure internal organs.

Do not modify the system. Use only the tools provided with the set and adhere to the instructions for use.

#### Risk of suffocation

This product contains small parts.

Keep small parts away from people who might swallow them, for example, small children.

#### Storage

Keep the product dry and away from sunlight. Store the product between 5 and 45 °C.

#### WARNING Risk of hyperglycaemia and infection

The use of sterile components after their Use by date can lead to infections and to abscesses.

Do not use components that have passed their Use by date.

#### Before You Get Started

Do not use this product for any medication other than U100 short-acting regular human insulin or sterile physiological saline solution.

Flush the new catheter with sterile physiological saline solution before use to make sure its tip is free.

#### It may be advisable to take a 30 min warm bath prior to the procedure to dissolve leftovers or crystallisation from insulin.

Distributed in the United Kingdom by: Roche Diabetes Care Limited, Burgess Hill, West Sussex, RH15 9RY, United Kingdom

Accu-Chek Pump Careline 1 UK Freephone number: 0800 731 22 91

1 calls may be recorded for training purposes

Some mobile operators may charge for calls to these numbers.

burgesshill.insulinpumps@roche.com

www.accu-check.co.uk

Australia Roche Diabetes Care Australia Pty. Limited

Insulin Pump Support: 1800 633 457

australia.insulinpumps@roche.com

www.accu-check.com.au

#### WARNING Risk of hypoglycaemia

When flushing the catheter or when inserting the guide wire during catheter replacement, the content of the catheter is delivered to the patient's body. The amount of insulin depends on the length of the implanted catheter:

- 9 cm (3.5 inches) catheter: 3.5 U (U100)
- 15 cm (5.9 inches) catheter: 6 U (U100)
- 25 cm (9.8 inches) catheter: 10 U (U100)

Closely monitor your patient's blood glucose during and after this procedure.

Introduce the guide wire into the catheter and push it through according to the catheter length.

4. Starting with the Seldinger technique, pull the catheter over the guide wire to remove it.

5. Select a new catheter with the appropriate length. Introduce the new catheter tip over the guide wire and move the catheter through the port body.

#### WARNING Risk of hyperglycaemia

A dislodged catheter hub can lead to insufficient insulin delivery. Make sure that the catheter hub is correctly aligned during and after the change of the membrane. The notch of the catheter hub fits into the port body in a defined direction only.

To position the catheter correctly, it is important to:

- Feed the catheter over the guide wire with the bend in the catheter opening towards the original implantation incision.
- Position the catheter to make sure that it fits correctly into the port body. Make sure that the catheter hub is correctly aligned before you insert it into the port body. The notch of the catheter hub fits into the port body in a defined direction only.

7. Use only the screw-in screwdriver with a pre-mounted membrane to screw the membrane into the port body. Screw in the pre-mounted membrane using a 90° angle and in one continuous, clockwise motion until it is hand-tight.

8. Use the disc applicator to position a new fixation disc on the port.

#### PRECAUTION Risk of hyperglycaemia

Handling errors or defective components can lead to insufficient insulin delivery.

After replacing a component of the system, check the blood glucose level at least once within 1 to 3 hours.

#### Care After Replacement

Stick the catheter replacement set sticker on the implant card carried by the patient.

Before commencing CIPI therapy, the catheter and the infusion set must be filled with the appropriate amount of insulin according to their lengths. If the patient's blood glucose level should unexpectedly increase, check for occlusions or leaks in the infusion set and review the log file of the insulin pump to make sure the catheter and the infusion set were filled according to the instructions for use.

#### WARNING Risk of hypoglycaemia

When filling the infusion set with insulin, the insulin pump delivers a considerable insulin amount on once.

Do not fill the infusion set while it is connected to the body. You may deliver an uncontrolled insulin amount.

#### PRECAUTION Risk of hyperglycaemia

The use of system components beyond the specified replacement intervals can cause occlusions or leaks that can lead to insufficient insulin delivery. Replace system components according to the specified intervals.

For more information, see the Accu-Chek DiaPort Handbook for Port Users.

#### Disposal

This set contains pointed or sharp-edged components. Dispose of pointed or sharp-edged components according to local regulations, in such a way that no third party can be injured. For more information, contact your local council or authority.

#### Reporting of Serious Incidents

For a patient/user/third party in the European Union and in countries with identical regulatory regime; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

#### Summary of Safety and Clinical Performance

You find the SSCP after the launch of the European database on medical devices (Eudamed) here: <https://ec.europa.eu/tools/eudamed>

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